® *Cyclotest Lady*

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Instruction manual

Mode d'emploi

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Natural Family
Planning

Planning familial
naturel

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Natural Family Planning

®Cyclotest Lady:

Special digital thermometer for measuring awaking temperature.

®Cyclotest Lady in daily use

_____ *battery compartment*

_____ *on/off button*

_____ *4 digit display*

_____ *flexible tip*

_____ *sensor*

When should you measure:

always in the morning **before** getting up –
at the same time if possible.

Where to take a measurement:

oral (under the tongue to the left or right of the
fraenum), vaginal or rectal – don't alterate during
a period.

When do you switch it off:

when a successful measurement has been acous-
tically indicated. (Forgotten to switch off? No
problem: ®Cyclotest Lady switches itself off auto-
matically after 10 minutes).

When do you make a record:

after getting up. The result stays in the memory
until the next measurement.

What to do with the result:

enter the measured value in the ®Cyclotest
period chart.

How to evaluate the results:

follow the instructions to be found on the follow-
ing pages.

You may congratulate yourself!

Like many other couples, you have decided to use the natural method of family planning and have acquired the °Cyclotest Lady precision device for this purpose.

°Cyclotest Lady is the ideal special thermometer for natural family planning. Even when this method demands a little more time and self-discipline – don't get discouraged and make it into a partnership! It has been shown to be a good thing when the partner actively takes part by entering the measured values into the period chart.

The fine variations in basal or awakening temperature give you an idea what is happening during your monthly cycle.

Normally only a single egg cell (ovum) which is capable of being fertilized is formed between two regular periods. Ovulation which occurs approximately in the middle of the cycle is decisive for both family planning and contraception. Medical research has shown that 1-2 days after ovulation, the basal temperature of a woman rises by 0.2 – 0.5°C. (Directly before the day of ovulation a drop of 0.1°C is regularly observed). The increased basal temperature stays at approximately the same level until the next period and then sinks again (see sample curve). This sequence is the normal one: significant deviations are a sign of irregularities of the period or pregnancy.

Use in pregnancy planning (family planning)

The unfertilized female egg cell (ovum) is only capable of being fertilized for several hours after ovulation, male sperm on the other hand is germinative for 2-3 days. For this reason, fertilization of the egg (conception) can only take place in the last 5 days before the above-mentioned temperature rise. If you can exactly determine the degree of fluctuation of your ovulation days over several cycles, then you have a good idea when you can expect your conception receptive days. Should the

GB-18

temperature stay high for more than 18 days after the initial rise, the expectations of a pregnancy are good.

Use in contraception

The fact that the female egg is only continuously viable for a few hours means that even the day after ovulation it will probably not let itself be fertilized. However, we naturally want to be sure that the recorded temperature rise is actually a result of ovulation and not caused by other factors such as stress, a cold, change of climate, lack of sleep etc.

According to Prof. Döring, one recognizes the „right“ temperature rise by the fact that the increase occurs within 48 hours or less and the temperature in the following three days is at least 0.2°C higher as in the previous 6 days. The first of the three days with increased temperature is then the temperature rise. This is followed by two periods of conception free days which will be explained in the following „temperature method“ paragraph.

a) The temperature method

This method is regularly confused with the „calendar method“ according to Prof. Knaus-Ogino which is no longer considered to be reliable. In the latter, the day of ovulation is more presumed than measured whereas in the temperature method it is measured. By using the special „[®]Cyclotest Lady“ thermometer, you can measure the fine variations of the basal temperature, record the value in the period chart and thus see when your possible conception and conception free periods are. With the „strict form“ of the temperature method (sexual intercourse only from the third day of the temperature rise till the onset of menstruation), you achieve almost the same reliability as the pill – and this in a perfectly natural way without taking medication or using mechanical methods. As Prof. Döring in „Obstetrics and Female Therapeutics“ 4/85 points out, the „strict form“ of the temperature method (Pearl-Index 0.8) and the pill (Pearl-Index 0.25-0.5) are right at the top of the reliability scale.

Also the „extended form“ of the temperature method is still as reliable as mechanical contraception methods. Here you can also have sexual intercourse in the time from the begin of the period up to six days before the temperature rise. To be precise: calculate 6 days back from the earliest temperature rise which you have determined from at least 6 sequential cycles.

b.) The symptothermal method

In this conception, „thermal“ stands for the previously described temperature method and is combined here with the observation of a „symptom“, of secretion at the neck of the uterus. In particular, menstruation is usually followed by a phase of infertile „dry days“ during which no so-called cervical secretion can be observed. A few days before the fertile phase, the secretion causes a moist and sticky feeling at the entrance to the vagina, the characteristics of which are thread-likeness, transparency and strong liquefaction. If you can make this observation 3-4 days before a (possibly only flat) temperature rise, then you can rely on the fact that this temperature rise is the correct one, ovulation has occurred.

When using the „[®]Cyclotest Lady“, please observe the following points:

■ Start records on the first day of menstruation – use a new period chart. The difference between period and irregular bleeding is a typical fall in temperature.

■ Carry out measurements daily: in the morning before getting up, if possible at approximately the same time. Don't eat anything before taking the measurement and avoid physical activity. In the case of inadequate sleep (at least 5-6 hours) omit the measurement and enter no value. As soon as you have determined the degree of fluctuation of the onset of your ovulation, you have enough information for family planning or contraception purposes if you take measurements every month starting shortly before the earliest temperature rise and continuing until the rise takes place as well as carrying out checks during menstruation.

■ Always measure in the same place: either in the mouth, in the vagina or the anus. Don't change within one cycle!

■ Enter the measured values into the period chart. The decisive indication of ovulation is, however, always the typical temperature rise.

■ Make a note of anything special: Make a note on the day concerned of events which could affect the temperature – intercourse (I), feverish cold (F), taking of medicine (M) or sleeping tablets (S), too little sleep (IS), change of climate (C – important!), alcohol (A), strain (S) and physical activity (P).

■ Also record bleeding (and condition of secretion): regular bleeding (P), strong, moderate, little irregular bleeding (sB, mB, lB), trace of blood (T); if possible also the occurrence and type of cervical secretion.

■ In the case of unusual events, give your doctor the period charts. We particularly recommend this should irregularities of your period occur. In case of doubt you should consult your doctor.


® **Cyclotest period chart (sample)**

The entries in a period chart are significant for the determination of the fertile phase. Carefully filled-in period charts supply important diagnostic indications and should always be available when visiting the doctor.

The counting of the cycle days commences with the onset of menstruation. Don't forget the month and calendar day.

The awakening or basaltemperature is entered daily.

Under (Symptoms) you can enter special observations.

Instructions for use

®Cyclotest Lady display

First place thermometer. Press the on/off button. You will here a short beep, then the self-check appears

followed by the last stored measured value indicated by the "M" symbol in the bottom right corner.

Subsequently the device indicates an internal test value of 37.00°C ($\pm 0.02^\circ\text{C}$). After this the device is ready for measurement.

A flashing „C" shows that a new measurement has started. As soon as the sensor has recognized a maximum value, the ®Cyclotest Lady emits a beep tone and the „C" ceases to blink. The measured value remains stored until the next reading.

®Cyclotest Lady – how does it work?

As a maximum thermometer, ®Cyclotest Lady follows the temperature rise (from room to body temperature) and retains the maximum value reached in its memory.

This procedure needs a maximum of 30 sec. After running for 10 minutes the ®Cyclotest Lady switches itself off automatically in order to save the battery.

Under normal use the battery has a life of 2 years. When the following character blinks, the battery has to be replaced.

This can be done by:

- Pulling the cap off
- Removing the battery cover and battery (possible with the help of a sharp object)

- Inserting new battery ensuring that the + sign is visible.
- Replace battery cover and cap.

Spare battery:

Silver-dioxide button cell Type L 41 or SR41,
1,55 V (Obtainable in shops selling watches)

Caution:

Batteries and technical devices do not belong in the domestic waste and must be deposited in the respective collection disposal point.

The measuring range:

®Cyclotest Lady can display temperatures from 32°C to 42.99°C. Should a value lie below or above this range, then „Lo °C“ or „Hi °C“ will be displayed respectively.

®Cyclotest Lady as fever thermometer

You can also use the ®Cyclotest Lady for measuring feverish temperatures. Look for the following characters in the display

36.59 °C

normal temperature

37.50 °C

high temperature or fever

Measurement Methods

Oral

Measuring time approx. 1 minute.

Vaginal/Rectal

Measuring time approx. 30 seconds.

Axillary

From the medical point of view this method always provides inaccurate values. And so it is unsuitable for accurate measuring.

Measuring time approx. 5-10 minutes.

Important Information

- Do not drop the device. It is not shock-proof.
- Before measurement check intactness of tip and display.
- This thermometer is determined only for measuring human body temperature.
- This thermometer contains small parts which can be swallowed by children. Don't leave children with this device uncontrolled.
- Do not leave the device in direct sunlight or near heat sources.
- Do not open the device (except the battery cover).
- For cleaning use only water or special substances recommended in chapter „Cleaning and Disinfection“.
- Cleansing bath max. 30 min.

Technical data

Type:	Maximum thermometer
Measuring range:	32.00 °C to 42.99 °C Temperature under 32.00 °C display "Lo" for low (too low) Temperature over 43.00 °C display "Hi" for high (too high)
Accuracy:	± 0.10 °C between 35.50 °C and 42.00 °C with surrounding temperature between 10° C and 40 °C relative air humidity 30-85 %
Display:	Liquid crystal display (LCD) with 4 digits, smallest display unit 0.01 °C
Acoustic Signal:	when temperature increase less than 0.1 °C within 16 sec
Storage temperature:	-10 °C to +60 °C relative air humidity 25-90 %
Battery:	Type LR41 / SR41, 1.55 V
Battery check:	Optic display ("  ")

Degree of protection against electric shock:
Type BF

Observe the operating instructions

This thermometer conforms to the
Directive for Medical Devices (93/42/EWG)
DIN EN 12470-3: 2000 Compact electrical ther-
mometers with maximum temperature indicator.

Manufacturer's guarantee

The equipment was manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we are prepared to give a guarantee in accordance with the following conditions:

During the guarantee period of 2 years from the date of purchase, we reserve the right to either repair any such defect at our expense or supply a perfect replacement unit (in any case the defective equipment must be returned to our factory).

Excluded from the guarantee are parts subject to normal wear and tear as well as damage caused by non-compliance with the operating instructions, improper handling (for example: unsuitable supply mains, breakages, leaking batteries) and/or dismantling of the unit by the purchaser. Furthermore, no claims for damages against us can be justified through the guarantee.

Claims under the guarantee can only be made during the guarantee period by presentation of the purchase receipt. In the case of claim under the guarantee, the equipment together with the purchase receipt and a written justification of the claim are to be sent to

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In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the supplier/salesman in accordance with § 437 BGB (The German Civil Code) are not limited by the guarantee.

**Conception?
Contraception?
®Cyclotest!**

®Cyclotest diagnosis products
have become indispensable for
natural family planning (NFP).

®Cyclotest 2 Plus

Cycle temperature monitor
Your path to your planned child

®Cyclotest BASIC

Cycle temperature monitor
Contraception without hormones

®Cyclotest

Special thermometer
with straddled scale for indication of mini-
mum temperature differences

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  Copyright 2006 UEBE Medical GmbH
Warenzeichen der UEBE Medical GmbH